

Stages and Platforms

Buyer's Guide

Scientifica's wide range of stages and platforms offer a variety of options for mounting your microscope, sample and other equipment.

They ensure superb stability and smooth movement for the most demanding electrophysiology experiments or advanced imaging applications



Discover more

MSP




UMS



MMTP



	Suitable Microscope	Movable	Motorised	In Vitro	In Vivo	Travel	Advanced Imaging
 <u>Motorised Movable Base Plate (MMBP)</u>	Upright	Yes	Yes	Yes	Yes	X Travel: 50mm Y Travel: 50mm	Yes
 <u>Motorised Movable Top Plate (MMTP)</u>	Upright	Yes	Yes	Yes	No	X Travel: 25mm Y Travel: 25mm	Yes
 <u>Movable Top Plate (MTP)</u>	Upright	Yes	No	Yes	No	X Travel: 25mm Y Travel: 25mm	Yes
 <u>SlicePlatform</u>	Upright	No	No	Yes	No	N/A	No
 <u>Motorised Movable Sample Plate (MMSP)</u>	Upright	Yes	Yes	Yes	Yes	X Travel: 20mm Y Travel: 20mm	Yes
 <u>Movable Sample Plate (MSP)</u>	Upright	Yes	No	Yes	Yes	X Travel: 20mm Y Travel: 20mm	Yes
 <u>Sample Post & Platform</u>	Upright	No	No	Yes	Yes	N/A	No

	Suitable Microscope	Movable	Motorised	In Vitro	In Vivo	Travel	Advanced Imaging
	Inverted	Yes	Yes	Yes	No	X Travel: 25mm Y Travel: 25mm	Yes
<u>Inverted Motorised Movable Top Plate (IMMTP)</u>							
	Inverted	Yes	No	Yes	No	X Travel: 25mm Y Travel: 25mm	Yes
<u>Inverted Movable Top Plate (IMTP)</u>							
	Inverted	No	No	Yes	No	N/A	No
<u>Inverted Fixed Platform</u>							
	Upright/ Inverted	No	No	Yes	Yes	N/A	No
<u>Manipulator Post and Platform</u>							
	Upright/ Inverted	Yes	No	Yes	No	X Travel: 25mm Y Travel: 25mm	No
<u>Manual XY Stage</u>							
	Upright/ Inverted	Yes	Yes	Yes	No	X Travel: 50mm Y Travel: 50mm	No
<u>Universal Motorised Stage (UMS)</u>							

Discover
more

